

# Infrared Analysis of Thin Multilayered Polymer Film Using Cantilever Enhanced Photoacoustic Detector

Kaori Watanabe and Fumiko Hiwatashi 

Systems Engineering Inc. [www.systems-eng.co.jp](http://www.systems-eng.co.jp) [info@systems-eng.co.jp](mailto:info@systems-eng.co.jp)

## Introduction

### Polymer thin film is used in...

- Structural material
- Electronic device
- Adhesive
- Paint and printing
- Biomaterial

nm to  $\mu\text{m}$  ordered polymer structure  
Multilayered films for hybrid and novel functions  
**Depth profiling analysis is necessary**


### Surface analysis

- XPS, SSIMS, SPM, XRD etc.
- Vibrational spectroscopy -> definitive molecular information
  - » PAS: Non-destructive (little sample preparation), depth-profiling

### Novel PAS detector (GASERA; booth 555)



- Cantilever and laser interferometer based detector
- **Ten times more sensitive than traditional acoustic microphone**




### This study

- Infrared PAS analysis of multilayered films
- Determination of accurate sampling depth of PAS with cantilever detector

## Photoacoustic Spectroscopy (PAS)

### Photoacoustic spectroscopy


1. Sample absorbs incident light
2. Absorbed energy partially transformed into heat (non-radiating relaxation)
3. Heat diffuses back to sample surface
4. Thermal expansion of surrounding gas detected



### PAS advantage

- Highly sensitive
- Non-destructive
- Depth-profiling
- Versatile
  - » Irregular shaped, small, or opaque sample etc.
  - » THz, IR, NIR, UV-Vis etc.

## PAS Sampling Depth




$\mu$ : "PAS Sampling Depth"  
Generated heat decays to  $1/e$

$\alpha$ : Thermal diffusivity  $\alpha = k/\rho C$  [ $\text{cm}^2 \text{s}^{-1}$ ]     $\omega$ : Angular modulation frequency  $\omega = 2\pi Vv$   
 $k$ : Thermal conductivity [ $\text{W cm}^{-1} \text{K}^{-1}$ ]     $V$ : Mirror velocity [ $\text{cm s}^{-1}$ ]  
 $\rho$ : Density [ $\text{kg cm}^{-3}$ ]     $v$ : Infrared wavenumber [ $\text{cm}^{-1}$ ]  
 $C$ : Specific heat [ $\text{J kg}^{-1} \text{K}^{-1}$ ]

## Results and Discussion

### 1. PAS spectrum of multilayered film

Layered sheets of food wrap  
PE, PVC, PVDC; 10  $\mu\text{m}$  each  
Mirror velocity  $V = 5 \text{ kHz} = 0.32 \text{ cm s}^{-1}$




1740  $\text{cm}^{-1}$  peak of PVC observed  
Sampling depth  $> 10 \mu\text{m}$   $\leftrightarrow$  PAS sampling depth  $\mu = 7.6 \mu\text{m}$

Photoacoustic signal observed beyond  $\mu$

### 2. PAS spectrum of multilayered film, various mirror velocity

Multilayered polyimide film  
Kapton® 300F929  
PI 50  $\mu\text{m}$  Teflon® FEP 12.5  $\mu\text{m}$




### Estimation of accurate sampling depth $A\mu$ (A: factor)

- 1775  $\text{cm}^{-1}$  peak: Observed at 10 kHz, not observed at 20 kHz
  - »  $3.8A < 12.5 \mu\text{m} < 5.3A$
- 1080  $\text{cm}^{-1}$  peak: Observed at 10 kHz, not observed at 20 kHz
  - »  $3.4A < 12.5 \mu\text{m} < 4.8A$
- 605  $\text{cm}^{-1}$  peak: Observed at 20 kHz, not observed at 40 kHz
  - »  $4.6A < 12.5 \mu\text{m} < 6.4A$

Estimated factor

$$2.6 < A < 2.7$$

i.e.,  $A \approx e$



Heat generated at depth  $\epsilon\mu$  can be detected  
Generated heat decays to  $1/e^2$  at surface

## Conclusion

Photoacoustic signal can be detected beyond PAS sampling depth  $\mu$ .  
Accurate sampling depth is estimated to be  $\sim \epsilon\mu$ , where the generated signal decays to  $1/e^2$  upon reaching the sample surface.